



LEG / ANKLE / FOOT POSITIONING



Heel Support

Reusable polymer heel support designed to float heels above table pad surface while in supine position.

Soft heel support features a v-shaped notch to position and accommodate the Achilles tendon and sculpted heel cup to cradle and securely position the foot in a neutral position.

40502 6 L x 3 1/4 W x 1 3/4" 15 x 8.3 x 4.5 cm



Heel Protector Pad

This "clamshell" style pad provides full coverage ankle and heel protection during lengthy surgeries. Interior surface of the pad is permanently labeled to ensure proper foot placement. Easy to use and secure with adjustable hook and loop attachment.

Accommodates approximately, a 15" ankle circumference. Use alone in supine and lateral position, or combine with Ankle Strap Stirrup and Candy Cane pads for lithotomy position.

40401 12 1/2 L x 7 3/4 W x 1/2 H" 32 x 20 x 1.2 cm



Dome-Shaped Positioners

Semi-circular design with flat bottom configuration provides pressure protection under the knees while the patient is in a supine position.

Helps to protect patient against foot drop for patients in prone position.

Can be used under the metatarsal or under the foot, against the board, or in the semi-Fowler position.

40603 14 L x 6 W x 3 H" 36 x 15 x 7.6 cm

40603L 20 L x 6 W x 3 H" 50.8 x 15 x 7.6 cm



Shown with Heel Support and Overlay

Knee Crutch Pads

Three sizes and thicknesses available to line knee crutches and help reduce pressure and shear in the popliteal space.

SMALL: Thinnest pad, best for lightweight support for urology, OB/ GYN, and pediatric procedures.

MEDIUM: 85% larger than the Small, this 1/4" pad offers midweight support for urology and OB/GYN procedures.

LARGE: Heavyweight 1/2" thick, 22 x 20" pad best suited for large knee crutches.

40402 Small 11 L x 9 W x 3/8 H" 27.9 x 22.9 x .95 cm

40403 Medium 17 L x 11 W x 1/4 H" 43 x 27.9 x .6 cm

40404 Large 22 L x 20 W x 1/2 H" 56 x 50 x 1.2 cm



Action Products, Inc.

800.228.7763 • 301.797.1414 • Fax 877.733.2073
www.actionproducts.com • service@actionproducts.com

ACTION and AKTON are registered trademarks of Action Products, Inc. © 2024 Action Products, Inc.

ACTION products feature our proprietary Akton® polymer



Toboggan Pad

This pad is for use with the arm/ leg guard shield (toboggan). Together they secure the arm or leg position on the table and protect the patient from being leaned on or pinched during the procedure.

The polymer pad cushions and helps protect the arm's ulnar and radial nerves.

The small film pocket slips on to the shield and pad is secured by hook & loop attachment (included).

The shield is only sold with the set.



Shown with Heel Support



Hook and loop attachment



40130	Toboggan Pad	16 ¼ L x 8 ¼ W x ¾ H"	41.27 x 20.96 x .95 cm
S40130*	Toboggan Pad Set	16 ¼ L x 8 ¼ W x ¾ H"	41.27 x 20.96 x .95 cm
	Shield size =	15 ½ L x 18 ¼ W"	39.37 x 46.4 cm

*Warranty covers pad only.

*Toboggan Armguard Shield is not available for sale or shipment outside of the United States and Canada.



Candy Cane Stirrup Pad Sets


Use in lithotomy position during gynecological and urological procedures; protects peroneal nerve from pressure and shear injuries that can result from leg contacting metal stirrups. Soft, durable and reusable product eliminates need for tedious wrapping and rewrapping of the leg holder posts. Recommend using Action® Ankle Stirrup Pad Set and full coverage Heel/ Ankle Protectors for enhanced patient protection.

**Lithotomy holders not included*

Ankle Strap Stirrup Pad Sets

Two adjustable polymer straps are combined to support and protect the patient's leg and foot while leg is suspended in lithotomy position.

40501	Candy Cane (Set of 2)	28 L x 2 OD x 1 ID"	71 x 5 x 2.5 cm
40510	Ankle Strap Stirrup Pads, Small (Set of 2)	12 ½ L x 1 ¾ W x ¾ H"	31.7 x 4.7 x .95 cm
40511	Ankle Strap Stirrup Pads, Large (Set of 2)	17 ¾ L x 1 ¾ W x ¾ H"	44.1 x 4.7 x .95 cm



Heel Comfort Cup Set

Lightweight polymer heel cups for comfort and support protects heels from pressure and shear injury. Can be used in conjunction with anti-embolism stockings. Appropriate for recovery and radiology. Comfortable all-polymer fastening strap attaches with simple hook and loop. One size fits all.

40520	(Set of 2)	4 ¼ L x 3 ¾ W x ½ H"	10.8 x 9.2 x 1.2 cm
--------------	------------	----------------------	---------------------



Heel & Elbow Protectors

Recommended for those clients "at risk" or diagnosed with compromised skin disorders. Full polymer strap offers enhanced, adjustable protection for heel or elbow.

20401	1/2" Thick
20401A	1/4" Thick

Positioning Strap

Welded seams and fluid-resistant film combined with a soft patient-facing polymer surface make this a must have for the operating room.

Radiolucent and MRI safe.

Hook and loop attachments secure the strap to standard tables. The longer strap was specifically designed for larger patients.

This strap is not meant to restrain patient.



40713	57 L x 4 W x ¼" H	144.78 x 10.16 x .64 cm
	AKTON polymer surface measures 98 x 8.2 x 0.6cm	
	Adjustable range: 28 x 43" 71.12 x 109.22 cm	
40713L	90 L x 4 W x ¼" H	228.6 x 10.16 x .64 cm
	AKTON polymer surface measures 190.5 x 8.2 x 0.6cm	
	Adjustable range: 59½ x 79½" 151.13 x 201.93 cm	



One Piece



Two Piece Set



Frog Leg Positioner

Designed for saphenous vein harvesting; the contouring and construction of the pad enhance site access while providing stable support.


Akton® polymer overlay protects from pressure and shear injury and stable foam base (visco foam atop high-density foam) is encased in 4-way stretch water resistant material that is fully welded to prevent fluid ingress.

Reusable product eliminates need for wasteful use of stacked blankets or towels and sandbags.

Replacement polymer overlay is available, Item 40718LPA. Overlay features hook and loop to attach it to the fabric-covered base.

Two-piece design allows for independent leg positioning.

40718LP	Frog Leg Positioner	32 L x 21 W x 5 H" includes 1/4" AKTON polymer	81.3 x 53.3 x 12.7 cm includes 0.6 cm AKTON polymer
40718LPA	Frog Leg Polymer Replacement		
40718LP-2	Frog Leg Positioner (2-Piece Positioner)	32 L x 10 ½ W x 5 H" includes 1/4" AKTON polymer	81.3 x 26.6 x 12.7 cm includes 0.6 cm AKTON polymer



Boot Stirrup Pad Set

Soft, reusable 1/4" thick polymer lithotomy boot liner designed to protect patient's lower leg and heel area by reducing pressure and shear.



Shown with SAM Boot

For use with universal stirrups and includes hook and loop attachments to secure to hardware. Used in lithotomy position for gynecological and urological cases.

40500	Universal (Set of 2)	23 ¾ L x 5 ⅝ W x ¼ H"	60.3 x 13.7 x .6 cm
40503	SAM (Set of 2)	27 L x 14 ¼ W x ¼ H"	68.5 x 36 x .6 cm

Fracture Table Foot Pad

For use on leg-prep stand and fracture table foot holders, for protection from pressure and shear injury.



40400	12 L x 4 W x ¼ H"	30.5 x 10 x .6 cm
--------------	-------------------	-------------------